



PG OTI-ONE

Electrical surgical table solutions



Stainless steel structure with slow rebound latex foam mattress
Five section body layout for multiple surgical position request
Reliable 6-movements electrical controlling system with smooth movement and silent lifting
Table top sliding for X-Ray and C-Arm
Individual kidney bridge elevation system
Built-in back up battery

OPERATING TABLE



PG OTI-ONE

PG OTI-ONE



TECHNICAL FEATURES



Designed specially to meet the different demands of general surgical operations, latex foam mattress, concealed rubber lift locker, single remote controller.

Class I, type B

Length: 2070mm

Width: 550mm

Height adjustment: 700-1000mm

Reversed Trendelenburg: 25°

Trendelenburg: 25°

Tilt left/right: +15°/15°

Horizontal sliding: 350mm

Net weight overall: 270Kg

Loading capacity: 250 kg

Leg plate up/down: +15°/-90°

Head rest up/down: +45°/-90°

Back plate up/down: +75°/-20°

Split leg plate: 90° (each)

Kidney bridge elevation: 120mm

Backup battery: included, for 80 movements

Power supply: 100-240V / 50-60Hz

Standard accessories: mattress pad (1 pc), cable controller (1 pc), arm support (1 pair), shoulder support (1 pair) body support (1 pair), leg support (1 pair), anesthesia frame (1 set)

PG OTI-ONE Plus



TECHNICAL FEATURES



Polyurethane mattress with sand spray surface, pedal wheel locker, double remote controller, one button reset, head plate angle adjustment, leg plate electrical control, 4 medical wheels (2 with brakes)

Class I, type B

Length: 2030mm

Width: 600mm

Height adjustment: 700-1000mm

Reversed Trendelenburg: 25°

Trendelenburg: 25°

Tilt left/right: +25°/25°

Horizontal sliding: 350mm

Net weight overall: 320Kg

Loading capacity: 300kg

Max loading under sliding process: 200kg

Leg plate electrical adjustment: +15°/-90°

Head plate adjustment: +40°/-90°

Back plate up/down: 70°/25°

Split leg plate: 200°

Kidney bridge elevation: 120mm

Dimensions: head plate: 300mm, back plate: 380mm seat

plate: 540mm, leg plate: 580mm, base size: 103x65x20cm

Medical wheel diameter: 125mm

Backup battery: included for 180 movements

Power supply: 100-240V / 50-60Hz

Standard accessories: anesthesia frame with IV pole (1 pc), arm support (1 pair), shoulder support (1 pair), lateral support (1 pair), leg support (1 pair), foot plate (1 pair)



progetti
Medical Equipment Solutions

PROGETTI S.r.l.
Strada del Rondello, 5
10028 Trofarello (Torino) - Italy
Ph. +39 011 644 738
Fax +39 011 645 822
info@progettimedical.com
www.progettimedical.com



ISO 13485:2016 CERTIFICATION



ISO 9001
BUREAU VERITAS
Certification
Quality system ISO certified

Distributed by: